

Pouring SDS-PAGE Gels

	Resolving				Stacking	
	<u>5%</u>	<u>7.5%</u>	<u>10%</u>	<u>12.5%</u>	<u>15%</u>	<u>4%</u>
dH₂O	6.1 mL	5.5 mL	4.9 mL	4.2 mL	3.6 mL	3.8 mL
40% Acrylamide/Bis	1.25 mL	1.9 mL	2.5 mL	3.1 mL	3.8 mL	0.5 mL
1.5M Tris-HCl pH 8.8	2.5 mL	2.5 mL	2.5 mL	2.5 mL	2.5 mL	-
1.0M Tris-HCl pH 6.8	-	-	-	-	-	0.64 mL
10% SDS	100 µL	100 µL	100 µL	100 µL	100 µL	50 µL
10% APS	50 µL	50 µL	50 µL	50 µL	50 µL	25 µL
TEMED	5 µL	5 µL	5 µL	5 µL	5 µL	5 µL
Total Volume (2 gels)	10 mL	10 mL	10 mL	10 mL	10 mL	5 mL

